



West Rigid Needle Shields for Staked-Needle, Pre-Filled Syringes

We've Got It Covered!



West Pre-Filled Syringe Components



Syringe plungers/pistons



Rigid and soft needle shields




Tip caps



Plunger rods



West Rigid Needle Shields: Protection & Safety



- Needle shields with rigid covers are recommended closures for prefilled syringes protecting the hypodermic needle tips through de-shielding via an optimized design.
- Combines the sealing properties of soft rubber shields with the rigidity of a polypropylene (PP) shell safeguarding against needle tip damage and deformation
- Manufactured from elastomers especially designed for prefilled syringe applications
- Helps assure container closure integrity through two sealing areas
- Designed to dimensionally fit ISO 11040-4 glass syringes with staked needles

* 7025/65 is not steam sterilizable in assembled syringe format



Formulation

- ✓ Formulation 7028/55¹ designed for needle shields
- ✓ High resistance to ozone cracking² helping ensure container closure integrity

- High resistance to ozone cracking²
- Designed for compliance with applicable EP and USP compendia
- Reduced risk of coring and fragmentation³
- Less susceptible to frosting, visual discoloration and blooming than traditional elastomers



In-Process Excellence

- ✓ Tight dimensional tolerances achieved through precision injection molding providing confidence
- ✓ In-process automated camera inspection resulting in dimensional reliability for processing

- Dimensional specifications with tight tolerances
- Designed for easy distribution and orientation on syringe assembly equipment from different vendors
- Utilizes existing assembly machine and filling equipment with minor adjustments
- Suitable for polymer (e.g., COP) syringes
- Usable with needles with safety features




Extensively Studied with Industry Partners

- ✓ Extensively studied in collaboration with industry partners to characterize³ reaffirming performance

- Reduced risk of pop-off during sterilization or pressurization³
- Ability to remain intact under high separation force³
- Pull-off force enabling controlled and seamless shield removal³

1. Formulation Characteristics, Elemental Extractables and Nitrosamines Data Sheets supported by Compliance Bulletins available upon customer request
 2. Technical Document: The Effects of Ozone and Various Sterilization Techniques on Elastomeric Rubber Needle Shields and Tip Caps
 3. TR 2018/183 West 5/8" RNS Product Characterization and TR 2018/182 West 1/2" RNS Product Characterization

Available in Two Sizes & Two Rubber Formulations*



RNS™ 1/2" (13mm) – Design 4144

RNS™ 5/8" (16mm) – Design 4145

*7028/55 and 7025/65 rubber formulations

7028/55 Rubber Formulation Details



Parameters	7028/55 Rubber Formulation ¹
Base Polymer	Synthetic Polyisoprene - Proprietary Elastomer Blend
Dry Natural Latex Rubber	Not Used As Raw Material
EP 3.2.9 chem. tests, USP <381> physicochem. tests	Complies to Type I
CHP & YBB Compliance	Current China National Standards ²
Coring	Low ³
Gas & Vapor Permeability	High ³
Frosting & Blooming	Low ⁴
Sterilization Compatibility	Steam, EtO and Gamma Irradiation ^{5, 6}
Resistance to Ozone Cracking	High ⁶
Extractables & Leachables (E&L)	Low ⁴
Extractable Elements	Zinc ⁷

1. Formulation Characteristics, Elemental Extractables and Nitrosamines Data Sheets supported by Compliance Bulletins available upon customer request
2. Not applicable for all part codes; request for China letter of authorization process
3. TR 2018/183 West 5/8" RNS Product Characterization & TR 2018/182 West 1/2" RNS Product Characterization
4. Data on file at West
5. Steam sterilization is not recommended for assembled prefilled syringes with a staked needle due to the adhesive used to assemble staked needles; rigid cover has potential of color change after gamma irradiation
6. Technical Document: The Effects of Ozone and Various Sterilization Techniques on Elastomeric Rubber Needle Shields and Tip Caps
7. Alternative rubber formulation 7025/65 offered with lower zinc levels; reference the Elemental Extractables Data Sheet for 7028/55 and 7025/65

Difference Is In West Manufacturing Process



**The difference is in
the West
manufacturing
process**

- West has a process to manufacture (soft) needle shields using a Precision Injection Molding (PIM) Cell approach
 - Approach is a combination of PIM and primary vision inspection
 - PIM molding technology yields components with tight tolerances and precise geometry
 - Automated visual inspection offers industry leading quality

- Siliconization is tailored for the rubber formulation and design
 - Constant wash load process obtains a consistent silicone oil application within a target output range

- Frequent manufacturing in-process-quality controls (IPC) ensure exact components
 - IPC release criteria incl. dimensions, siliconization level, particulate, bioburden, and endotoxin

West Rigid Needle Shield Specifications



West **RNS**[™]
Rigid Needle Shields



- Visual
 - According to AQL calculation
- Pull Off Force
 - Pull-Off-Forces $5\text{N} \leq X \leq 35\text{N}$
- Standard Dimensions
- Endotoxin level:
 - $\leq 0.25 \text{ EU/cm}^2$
- Bioburden level:
 - $\leq 2.5 \text{ CFU/10 cm}^2$
- Particle Counting Index [PCI]
 - $\leq 6 / 10 \text{ cm}^2$

RNS™ Rigid Needle Shields Performance Evaluation



- Extensive study in cooperation with industry professionals
- Test Program: Functional testing over time; Machinability; Transport Simulation

Functional testing on all design-rubber combinations:

Pull-off Force [N]	Separation Force [N]	Coring	Pop-Off/Movement	Leak Test
<ul style="list-style-type: none"> • Important for proper de-shielding of the RNS from the syringe by patients or caregivers. • Pull-off-Force: 5-35 N 	<ul style="list-style-type: none"> • The rigid shell must be well joined with the flexible rubber component to guarantee a safe handling when disengaging the RNS from the syringe • Removal force ≥ 50 N 	<ul style="list-style-type: none"> • During RNS assembly on the syringe the needle may punch out fragments from the rubber component. Fragments represents a critical risk for the injection. 	<ul style="list-style-type: none"> • Simulation of steam sterilization and air transport which bares the risk of RNS movement on the syringe or pop off leading to the loss of integrity. • RNS Movement $\leq 0,5$ mm 	<ul style="list-style-type: none"> • To ensure integrity of assembled RNS syringes are filled with Wfl and an increasing pressure is applied to the plunger and the leakage is tracked.
Ethylene Oxide Sterilization	EtO Residual			
<ul style="list-style-type: none"> • Verify the EtO sterilizability of RNS assembled syringes with bioindicator spore strips. The rubber formulation and the windows of the rigid shell allow an appropriate gas/vapor penetration 	<ul style="list-style-type: none"> • Levels of remaining EtO and ethylene chlorodyne were tested compliant with ISO 10993-7 			


We've got it covered!



West **RNS**[™]
Rigid Needle Shields



Thank You!



This presentation is for informational purposes only. West's products and solutions are sold on the basis that it is the customer's responsibility to evaluate and test the West product or solution to determine its compatibility with other materials and fitness for any end use. WEST MAKES NO WARRANTIES, WHETHER EXPRESS, IMPLIED OR STATUTORY, INCLUDING, WITHOUT LIMITATION, THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE, RELATING TO THE INFORMATION IN THIS PRESENTATION.

West and the diamond logo and RNS are trademarks or registered trademarks of West Pharmaceutical Services, Inc. in the United States and other jurisdictions.

© 2024 by West Pharmaceutical Services, Inc., Exton, PA. All rights reserved.

This material is protected by copyright. No part of it may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying or otherwise, without written permission of West Pharmaceutical Services, Inc.