



By your side  
for a healthier world™



Leading the way with  
integrated containment  
and delivery solutions

## SMARTDOSE® 3.5 PRELOADED Injector

PRECISE. RELIABLE.  
PRELOADED.

### Expanded Platform

- SmartDose® Preloaded Delivery Platform
- Uses proprietary prefillable cartridge as primary container
- Final assembly does not require cleanroom environment

### Wearable Injector

- Reduced user steps with preloaded cartridge
- Subcutaneous injection for self-administration
- Pre-programmable injection time and rate
- Ability to deliver high volume and high viscosity drug products

### Patient-Centric

- Simple, easy-to-use system with fewer user steps
- Visual, tactile and audible feedback boosts patient confidence
- Onboarding and training solutions



-  **PRECISE**
-  **RELIABLE**
-  **PRELOADED**



### Technical Specifications

Target Tissue	Subcutaneous
Max Container Volume	Up to 3.7mL
Sterilization	Device sterilization not required Ready to Use (RU) sterile container system
Container Type	Custom cartridge with cannula/needle shield and Flurotec® piston
Container Material	Daikyo Crystal Zenith® cyclic olefin polymer (COP)
Tub & Nest Format	8x4 configuration for standard filling equipment
Needle Safety System	Automatic needle insertion and retraction
Driving Mechanism	Electromechanical
System Control	Pre-programmable
Cartridge Loading	Preloaded
Wireless Connectivity	Bluetooth
Temperature Sensing	Optional capability

For more information visit [www.westpharma.com/SmartDose](http://www.westpharma.com/SmartDose)

**North America** 1-800-345-9800 option 8 | **South America** +55 11 4055 6061 | **Europe** +49 (0) 2403-7960 | **Asia Pacific** +65 6862 3400

**West Pharmaceutical Services, Inc.** 530 Herman O. West Drive | Exton, PA 19341

West and the diamond logo and By your side for a healthier world are registered trademarks or trademarks of West Pharmaceutical Services, Inc., in the United States and other jurisdictions. SmartDose® is a registered trademark of West Pharma. Services IL, Ltd., a subsidiary of West Pharmaceutical Services, Inc. Daikyo Crystal Zenith® and Flurotec® are registered trademarks of Daikyo Seiko, Ltd. Daikyo Crystal Zenith® technology is licensed from Daikyo Seiko, Ltd.