



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

West Pharma. Services IL, Ltd.
(formerly Medimop Medical Projects)
West Pharmaceutical Services, Inc.

**2nd AMENDED URGENT: MEDICAL
DEVICE RECALL**

**VIAL2BAG (DC) 13mm
VIAL2BAG (DC) 20mm
VIAL2BAG
I.V. FLUID TRANSFER SETS**

February 01, 2019

To: West Pharma. Services Distributor / Customer,

West Pharma. Services IL, Ltd. (formerly Medimop Medical Projects) and West Pharmaceutical Services Inc. has expanded the voluntary recall for the product codes and lot numbers listed as follows:

Product Code	Product Description	Lot #
6070104	Vial2Bag	Reference Appendix #1
6070111	Vial2Bag Direct Connect (DC) 13mm	
6070112	Vial2Bag Direct Connect (DC) 20mm	
6070112C		

Reason for the Voluntary Recall:

West is recalling these products due to the possibility that the device may not adequately transfer the drug contents to the IV bag. If inadequate transfer occurs there is the potential that the device



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

may deliver variable or unpredictable dosing. **This notice should replace the URGENT AMENDED MEDICAL DEVICE RECALL dated January 24, 2018.**

We are aware of sixteen (16) products with complaints of a similar nature since March 8, 2017, reported with the use of oxytocin in conjunction with the Vial2Bag DC 13mm in pregnant women, in connection with labor and delivery.

Risk to Health:

The reports conveyed that healthcare professionals had observed effects of the drug oxytocin being unpredictable or variable when administered with the Vial2Bag DC 13mm. The events reported include severe and rapid uterine contractions (uterine tachysystole), deceleration in fetal heart rate (with subsequent cesarean section), and hemorrhage. West has not determined whether the observed variable or unpredictable dosing is limited to oxytocin delivery through the affected medical devices.

Variable and unpredictable dosing of medications, generally speaking, may lead to overdosing or under-dosing, with the effects depending upon the specific drug, the specific patient, and the condition for which the drug was ordered. Adverse health consequences may range from minor to serious injury and/or death.

There have been no reported events for any Vial2Bag, Vial2Bag DC 13mm, Vial2Bag DC 20mm product with an outcome of death.

Actions to be taken by the Customer/User:

Our records indicate that you have received products consisting of the Lots as listed in Appendix #1 (the "Affected Product") that are subject to this voluntary recall.

1. If you have any affected product in your inventory, please immediately discontinue use, remove it from your inventory, and quarantine the product to prevent inadvertent use.
2. Contact Progressive Medical Inc. (PMI), 997 Horan Drive Fenton, Missouri 63026-2401 to schedule return of Affected Products.
 - 1-800-969-6331; 8am to 5pm, Central Time, Monday through Friday
3. ***Please complete the Recall Acknowledgement and Receipt Form*** and return to Linda Bagby by email at cs@progressivemedinc.com when complete responses are available (See attached form). This will allow us to document your receipt of this letter.

4. Alternative Options to Admix Medication:

We recommend that you rely on the current standard of practice of medicine for admixing medications for I.V. therapy.

Actions to be taken by the Distributor

1. Please examine your inventory and quarantine product subject to recall.
2. If the product was distributed to your customers, please identify your customers.



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

3. Notify and provide this letter to all accounts/customers where the Affected Product was distributed.
4. Follow up with each account/customer to locate and identify all Affected Product per the attached recall return response document and establish the number of units of Affected Product in possession of each account/customer.
5. Please complete the **Recall Acknowledgement and Receipt Form** below and return (via email) to recall@westpharma.com when complete responses are available (See attached form).

REFUND / REPLACEMENT FOR RETURNED PRODUCT -

Returned product will be credited by the distributor.

OTHER INFORMATION:

You should report any adverse events or quality problems with these devices to the manufacturer, West Pharma. Services IL, Ltd. You can also report an adverse event or quality problem experienced with the use of this product to FDA's MedWatch Adverse Event Reporting program:

- Online at: <http://www.fda.gov/Safety/MedWatch/HowToReport/default.htm> or
- By telephone: 1-800-FDA-1088

West regrets any inconvenience that this action may cause and appreciates your understanding and cooperation. This action is being undertaken to ensure the highest level of patient safety and customer satisfaction.

Customers who have questions or need additional information regarding this recall may contact West's U.S. distributor Progressive Medial Inc. at 1-800-969-6331; 8am to 5pm, Central Time, Monday through Friday or at cs@progressivemedinc.com.

This letter can also be found on the West Pharmaceutical Services, Inc. website at: www.westpharma.com/support/alerts-and-notice

Sincerely,

DocuSigned by:
Andy Polywacz
9AEC5BD85E2F4B1...

01-Feb-2019 | 3:58 PM EST

Andrew Polywacz
Vice President Quality & RA



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

**2ND AMENDED URGENT MEDICAL DEVICE RECALL
ACKNOWLEDGEMENT AND RECEIPT FORM
Response is Required**

**Vial2Bag
Vial2Bag DC 13mm
Vial2Bag DC 20mm**

Product Number(s): 6070104, 6070111, 6070112, 6070112C

Check the appropriate box and return this form

I have read and understand the recall actions to be taken provided in the February 1, 2019 Letter.

We have no inventory within the scope of this recall.

We have the following Affected Product at our facility and have discontinued use and distribution. The Affected Product has been quarantined, and the following quantities shall be returned.

Product Code	Lot Number	Quantity (units)



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Please Print Legibly

(Print Name)	(Date)
(Signature)	(Telephone Number)
(Institution Name)	(Email Address)
(Institution Street Address)	Alternate Mailing Address
(Institution City, State, Zip)	(Street Address)
(Country)	(City, State, Zip)



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Appendix #1:

Vial2Bag			
Part#	Lot	UDI	Expiry Date
6070104	7885	NA	Mar-19
6070104	7889	NA	Mar-19
6070104	8300	(01)07290108240047(17)190601(10)8300(90)36098093(91)6070104	1-Jun-19
6070104	8391	(01)07290108240047(17)190701(10)8391(90)36098093(91)6070104	1-Jul-19
6070104	8646	(01)07290108240047(17)191001(10)8646(90)36098093(91)6070104	1-Oct-19
6070104	8788	(01)07290108240047(17)191201(10)8788(90)36098093(91)6070104	1-Dec-19
6070104	8950	(01)07290108240047(17)200101(10)8950	1-Jan-20
6070104	9038	(01)07290108240047(17)200301(10)9038(90)36098093(91)6070104	1-Mar-20
6070104	9039	(01)07290108240047(17)200401(10)9039(90)36098093(91)6070104	1-Apr-20
6070104	9166	(01)07290108240047(17)200401(10)9166(90)36098093(91)6070104	1-Apr-20
6070104	9167	(01)07290108240047(17)200401(10)9167(90)36098093(91)6070104	1-Apr-20
6070104	9669	(01)07290108240047(17)201001(10)9669(90)36098093(91)6070104	1-Oct-20
6070104	9682	(01)07290108240047(17)200901(10)9682	1-Sep-20
6070104	9982	(01)07290108240047(17)210601(10)9982(90)36098093(91)6070104	1-Jun-21
6070104	A672	(01)07290108240047(17)210601(10)A672(90)36098093(91)6070104	1-Dec-21
6070104	C070	(01)07290108240047(17)210601(10)C070	1-Dec-21

Vial2Bag DC 13mm			
Part#	Lot	UDI	Expiry Date
6070111	7688	(01)07290108240061(17)190101(10)7688(90)36098125(91)6070111	01-Jan-19
6070111	7752	(01)07290108240061(17)190201(10)7752(90)36098125(91)6070111	01-Feb-19
6070111	7909	(01)07290108240061(17)190301(10)7909(90)36098125(91)6070111	01-Mar-19
6070111	7955	(01)07290108240061(17)190501(10)7955(90)36098125(91)6070111	01-May-19
6070111	8152	(01)07290108240061(17)190801(10)8152(90)36098125(91)6070111	01-Aug-19
6070111	8153	(01)07290108240061(17)190801(10)8153(90)36098125(91)6070111	01-Aug-19
6070111	8154	(01)07290108240061(17)190801(10)8154(90)36098125(91)6070111	01-Aug-19
6070111	8232	(01)07290108240061(17)190901(10)8232(90)36098125(91)6070111	01-Sep-19
6070111	8233	(01)07290108240061(17)190901(10)8233(90)36098125(91)6070111	01-Sep-19
6070111	8234	(01)07290108240061(17)190901(10)8234(90)36098125(91)6070111	01-Sep-19
6070111	8235	(01)07290108240061(17)190901(10)8235(90)36098125(91)6070111	01-Sep-19
6070111	9028	(01)07290108240061(17)200301(10)9028(90)36098125(91)6070111	01-Mar-20
6070111	9029	(01)07290108240061(17)200301(10)9029(90)36098125(91)6070111	01-Mar-20
6070111	9175	(01)07290108240061(17)200401(10)9175(90)36098125(91)6070111	01-Apr-20
6070111	9176	(01)07290108240061(17)200501(10)9176(90)36098125(91)6070111	01-May-20
6070111	9177	(01)07290108240061(17)200501(10)9177(90)36098125(91)6070111	01-May-20
6070111	9178	(01)07290108240061(17)200501(10)9178(90)36098125(91)6070111	01-May-20
6070111	9367	(01)07290108240061(17)200701(10)9367(90)36098125(91)6070111	01-Jul-20
6070111	9368	(01)07290108240061(17)200701(10)9368(90)36098125(91)6070111	01-Jul-20



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 13mm			
Part#	Lot	UDI	Expiry Date
6070111	9369	(01)07290108240061(17)200701(10)9369(90)36098125(91)6070111	01-Jul-20
6070111	9370	(01)07290108240061(17)200701(10)9370(90)36098125(91)6070111	01-Jul-20
6070111	9655	(01)07290108240061(17)201001(10)9655(90)36098125(91)6070111	01-Oct-20
6070111	9656	(01)07290108240061(17)201001(10)9656(90)36098125(91)6070111	01-Oct-20
6070111	9710	(01)07290108240061(17)201101(10)9710(90)36098125(91)6070111	01-Nov-20
6070111	9711	(01)07290108240061(17)201101(10)9711(90)36098125(91)6070111	01-Nov-20
6070111	9712	(01)07290108240061(17)201201(10)9712(90)36098125(91)6070111	01-Dec-20
6070111	9713	(01)07290108240061(17)201201(10)9713(90)36098125(91)6070111	01-Dec-20
6070111	9839	(01)07290108240061(17)201201(10)9839(90)36098125(91)6070111	01-Dec-20
6070111	9840	(01)07290108240061(17)201201(10)9840(90)36098125(91)6070111	01-Dec-20
6070111	9918	(01)07290108240061(17)201101(10)9918(90)36098125(91)6070111	01-Nov-20
6070111	9983	(01)07290108240061(17)210101(10)9983(90)36098125(91)6070111	01-Jan-21
6070111	9984	(01)07290108240061(17)210201(10)9984(90)36098125(91)6070111	01-Feb-21
6070111	9992	(01)07290108240061(17)210501(10)9992(90)36098125(91)6070111	01-May-21
6070111	9993	(01)07290108240061(17)210601(10)9993(90)36098125(91)6070111	01-Jun-21
6070111	9994	(01)07290108240061(17)210601(10)9994(90)36098125(91)6070111	01-Jun-21
6070111	9995	(01)07290108240061(17)210601(10)9995(90)36098125(91)6070111	01-Jun-21
6070111	A123	(01)07290108240061(17)210701(10)A123(90)36098125(91)6070111	01-Jul-21
6070111	A124	(01)07290108240061(17)210701(10)A124(90)36098125(91)6070111	01-Jul-21
6070111	A790	(01)07290108240061(17)210801(10)A790(90)36098125(91)6070111	01-Aug-21
6070111	A791	(01)07290108240061(17)210901(10)A791(90)36098125(91)6070111	01-Sep-21
6070111	A792	(01)07290108240061(17)210901(10)A792(90)36098125(91)6070111	01-Sep-21
6070111	A793	(01)07290108240061(17)210901(10)A793(90)36098125(91)6070111	01-Sep-21
6070111	A794	(01)07290108240061(17)210901(10)A794(90)36098125(91)6070111	01-Sep-21
6070111	A795	(01)07290108240061(17)210901(10)A795(90)36098125(91)6070111	01-Sep-21
6070111	A796	(01)07290108240061(17)210901(10)A796(90)36098125(91)6070111	01-Sep-21
6070111	A797	(01)07290108240061(17)210901(10)A797(90)36098125(91)6070111	01-Sep-21
6070111	A805	(01)07290108240061(17)210801(10)A805(90)36098125(91)6070111	01-Aug-21
6070111	A806	(01)07290108240061(17)210801(10)A806(90)36098125(91)6070111	01-Aug-21
6070111	A807	(01)07290108240061(17)210901(10)A807(90)36098125(91)6070111	01-Sep-21
6070111	A808	(01)07290108240061(17)210901(10)A808(90)36098125(91)6070111	01-Sep-21
6070111	A809	(01)07290108240061(17)210901(10)A809(90)36098125(91)6070111	01-Sep-21
6070111	A810	(01)07290108240061(17)210901(10)A810(90)36098125(91)6070111	01-Sep-21
6070111	A834	(01)07290108240061(17)211001(10)A834(90)36098125(91)6070111	01-Oct-21
6070111	A835	(01)07290108240061(17)211001(10)A835(90)36098125(91)6070111	01-Oct-21
6070111	A836	(01)07290108240061(17)211001(10)A836(90)36098125(91)6070111	01-Oct-21
6070111	A837	(01)07290108240061(17)211001(10)A837(90)36098125(91)6070111	01-Oct-21
6070111	A838	(01)07290108240061(17)211001(10)A838(90)36098125(91)6070111	01-Oct-21
6070111	A839	(01)07290108240061(17)211001(10)A839(90)36098125(91)6070111	01-Oct-21



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	7865	(01)07290108240054(17)190201(10)7865(90)36098126(91)6070112	1-Feb-19
6070112	7866	(01)07290108240054(17)190201(10)7866(90)36098126(91)6070112	1-Feb-19
6070112	7867	(01)07290108240054(17)190201(10)7867(90)36098126(91)6070112	1-Feb-19
6070112	7868	(01)07290108240054(17)190201(10)7868(90)36098126(91)6070112	1-Feb-19
6070112	7886	(01)07290108240054(17)190301(10)7886(90)36098126(91)6070112	1-Mar-19
6070112	7887	(01)07290108240054(17)190301(10)7887(90)36098126(91)6070112	1-Mar-19
6070112	7888	(01)07290108240054(17)190301(10)7888(90)36098126(91)6070112	1-Mar-19
6070112	7900	(01)07290108240054(17)190301(10)7900(90)36098126(91)6070112	1-Mar-19
6070112	7901	(01)07290108240054(17)190401(10)7901(90)36098126(91)6070112	1-Apr-19
6070112	7902	(01)07290108240054(17)190401(10)7902(90)36098126(91)6070112	1-Apr-19
6070112	7949	(01)07290108240054(17)190501(10)7949(90)36098126(91)6070112	1-May-19
6070112	7951	(01)07290108240054(17)190501(10)7951(90)36098126(91)6070112	1-May-19
6070112	7952	(01)07290108240054(17)190501(10)7952(90)36098126(91)6070112	1-May-19
6070112	7954	(01)07290108240054(17)190501(10)7954(90)36098126(91)6070112	1-May-19
6070112	7957	(01)07290108240054(17)190401(10)7957(90)36098126(91)6070112	1-Apr-19
6070112	7958	(01)07290108240054(17)190401(10)7958(90)36098126(91)6070112	1-Apr-19
6070112	7959	(01)07290108240054(17)190401(10)7959(90)36098126(91)6070112	1-Apr-19
6070112	7960	(01)07290108240054(17)190401(10)7960(90)36098126(91)6070112	1-Apr-19
6070112	8019	(01)07290108240054(17)190601(10)8019(90)36098126(91)6070112	1-Jun-19
6070112	8020	(01)07290108240054(17)190601(10)8020(90)36098126(91)6070112	1-Jun-19
6070112	8021	(01)07290108240054(17)190601(10)8021(90)36098126(91)6070112	1-Jun-19
6070112	8022	(01)07290108240054(17)190601(10)8022(90)36098126(91)6070112	1-Jun-19
6070112	8132	(01)07290108240054(17)190601(10)8132(90)36098126(91)6070112	1-Jun-19
6070112	8133	(01)07290108240054(17)190601(10)8133(90)36098126(91)6070112	1-Jun-19
6070112	8134	(01)07290108240054(17)190601(10)8134(90)36098126(91)6070112	1-Jun-19
6070112	8135	(01)07290108240054(17)190701(10)8135(90)36098126(91)6070112	1-Jul-19
6070112	8136	(01)07290108240054(17)190701(10)8136(90)36098126(91)6070112	1-Jul-19
6070112	8137	(01)07290108240054(17)190701(10)8137(90)36098126(91)6070112	1-Jul-19
6070112	8138	(01)07290108240054(17)190701(10)8138(90)36098126(91)6070112	1-Jul-19
6070112	8139	(01)07290108240054(17)190701(10)8139(90)36098126(91)6070112	1-Jul-19
6070112	8140	(01)07290108240054(17)190701(10)8140(90)36098126(91)6070112	1-Jul-19
6070112	8141	(01)07290108240054(17)190701(10)8141(90)36098126(91)6070112	1-Jul-19
6070112	8142	(01)07290108240054(17)190701(10)8142(90)36098126(91)6070112	1-Jul-19
6070112	8143	(01)07290108240054(17)190701(10)8143(90)36098126(91)6070112	1-Jul-19
6070112	8144	(01)07290108240054(17)190701(10)8144(90)36098126(91)6070112	1-Jul-19
6070112	8145	(01)07290108240054(17)190701(10)8145(90)36098126(91)6070112	1-Jul-19



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	8147	(01)07290108240054(17)190701(10)8147(90)36098126(91)6070112	1-Jul-19
6070112	8148	(01)07290108240054(17)190801(10)8148(90)36098126(91)6070112	1-Aug-19
6070112	8149	(01)07290108240054(17)190801(10)8149(90)36098126(91)6070112	1-Aug-19
6070112	8150	(01)07290108240054(17)190801(10)8150(90)36098126(91)6070112	1-Aug-19
6070112	8151	(01)07290108240054(17)190801(10)8151(90)36098126(91)6070112	1-Aug-19
6070112	8470	(01)07290108240054(17)190801(10)8470(90)36098126(91)6070112	1-Aug-19
6070112	8471	(01)07290108240054(17)190901(10)8471(90)36098126(91)6070112	1-Sep-19
6070112	8472	(01)07290108240054(17)190901(10)8472(90)36098126(91)6070112	1-Sep-19
6070112	8473	(01)07290108240054(17)190901(10)8473(90)36098126(91)6070112	1-Sep-19
6070112	8483	(01)07290108240054(17)190901(10)8483(90)36098126(91)6070112	1-Sep-19
6070112	8553	(01)07290108240054(17)191001(10)8553(90)36098126(91)6070112	1-Oct-19
6070112	8554	(01)07290108240054(17)191001(10)8554(90)36098126(91)6070112	1-Oct-19
6070112	8555	(01)07290108240054(17)191001(10)8555(90)36098126(91)6070112	1-Oct-19
6070112	8556	(01)07290108240054(17)191001(10)8556(90)36098126(91)6070112	1-Oct-19
6070112	8557	(01)07290108240054(17)191001(10)8557(90)36098126(91)6070112	1-Oct-19
6070112	8558	(01)07290108240054(17)191001(10)8558(90)36098126(91)6070112	1-Oct-19
6070112	8559	(01)07290108240054(17)191101(10)8559(90)36098126(91)6070112	1-Nov-19
6070112	8560	(01)07290108240054(17)191101(10)8560(90)36098126(91)6070112	1-Nov-19
6070112	8561	(01)07290108240054(17)191101(10)8561(90)36098126(91)6070112	1-Nov-19
6070112	8562	(01)07290108240054(17)191101(10)8562(90)36098126(91)6070112	1-Nov-19
6070112	8563	(01)07290108240054(17)191101(10)8563(90)36098126(91)6070112	1-Nov-19
6070112	8564	(01)07290108240054(17)191101(10)8564(90)36098126(91)6070112	1-Nov-19
6070112	8565	(01)07290108240054(17)191101(10)8565(90)36098126(91)6070112	1-Nov-19
6070112	8566	(01)07290108240054(17)191101(10)8566(90)36098126(91)6070112	1-Nov-19
6070112	8567	(01)07290108240054(17)191101(10)8567(90)36098126(91)6070112	1-Nov-19
6070112	8568	(01)07290108240054(17)191101(10)8568(90)36098126(91)6070112	1-Nov-19
6070112	8569	(01)07290108240054(17)191201(10)8569(90)36098126(91)6070112	1-Dec-19
6070112	8570	(01)07290108240054(17)191201(10)8570(90)36098126(91)6070112	1-Dec-19
6070112	8571	(01)07290108240054(17)191201(10)8571(90)36098126(91)6070112	1-Dec-19
6070112	8696	(01)07290108240054(17)190901(10)8696(90)36098126(91)6070112	1-Sep-19
6070112	8760	(01)07290108240054(17)191201(10)8760(90)36098126(91)6070112	1-Dec-19
6070112	8761	(01)07290108240054(17)191201(10)8761(90)36098126(91)6070112	1-Dec-19
6070112	8762	(01)07290108240054(17)191201(10)8762(90)36098126(91)6070112	1-Dec-19
6070112	8763	(01)07290108240054(17)191201(10)8763(90)36098126(91)6070112	1-Dec-19
6070112	8764	(01)07290108240054(17)191201(10)8764(90)36098126(91)6070112	1-Dec-19
6070112	8765	(01)07290108240054(17)200101(10)8765(90)36098126(91)6070112	1-Jan-20
6070112	8766	(01)07290108240054(17)200101(10)8766(90)36098126(91)6070112	1-Jan-20



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	8767	(01)07290108240054(17)200101(10)8767(90)36098126(91)6070112	1-Jan-20
6070112	8768	(01)07290108240054(17)200101(10)8768(90)36098126(91)6070112	1-Jan-20
6070112	8769	(01)07290108240054(17)200101(10)8769(90)36098126(91)6070112	1-Jan-20
6070112	8770	(01)07290108240054(17)200101(10)8770(90)36098126(91)6070112	1-Jan-20
6070112	8771	(01)07290108240054(17)200101(10)8771(90)36098126(91)6070112	1-Jan-20
6070112	8772	(01)07290108240054(17)200101(10)8772(90)36098126(91)6070112	1-Jan-20
6070112	8887	(01)07290108240054(17)200101(10)8887(90)36098126(91)6070112	1-Jan-20
6070112	8888	(01)07290108240054(17)200201(10)8888(90)36098126(91)6070112	1-Feb-20
6070112	8889	(01)07290108240054(17)200201(10)8889(90)36098126(91)6070112	1-Feb-20
6070112	8890	(01)07290108240054(17)200201(10)8890(90)36098126(91)6070112	1-Feb-20
6070112	8891	(01)07290108240054(17)200201(10)8891(90)36098126(91)6070112	1-Feb-20
6070112	8892	(01)07290108240054(17)200201(10)8892(90)36098126(91)6070112	1-Feb-20
6070112	8893	(01)07290108240054(17)200201(10)8893(90)36098126(91)6070112	1-Feb-20
6070112	8894	(01)07290108240054(17)200201(10)8894(90)36098126(91)6070112	1-Feb-20
6070112	8895	(01)07290108240054(17)200201(10)8895(90)36098126(91)6070112	1-Feb-20
6070112	9040	(01)07290108240054(17)200301(10)9040(90)36098126(91)6070112	1-Mar-20
6070112	9041	(01)07290108240054(17)200301(10)9041(90)36098126(91)6070112	1-Mar-20
6070112	9042	(01)07290108240054(17)200301(10)9042(90)36098126(91)6070112	1-Mar-20
6070112	9043	(01)07290108240054(17)200401(10)9043(90)36098126(91)6070112	1-Apr-20
6070112	9044	(01)07290108240054(17)200401(10)9044(90)36098126(91)6070112	1-Apr-20
6070112	9045	(01)07290108240054(17)200401(10)9045(90)36098126(91)6070112	1-Apr-20
6070112	9046	(01)07290108240054(17)200201(10)9046(90)36098126(91)6070112	1-Feb-20
6070112	9047	(01)07290108240054(17)200201(10)9047(90)36098126(91)6070112	1-Feb-20
6070112	9169	(01)07290108240054(17)200401(10)9169(90)36098126(91)6070112	1-Apr-20
6070112	9170	(01)07290108240054(17)200501(10)9170(90)36098126(91)6070112	1-May-20
6070112	9171	(01)07290108240054(17)200501(10)9171(90)36098126(91)6070112	1-May-20
6070112	9172	(01)07290108240054(17)200501(10)9172(90)36098126(91)6070112	1-May-20
6070112	9173	(01)07290108240054(17)200501(10)9173(90)36098126(91)6070112	1-May-20
6070112	9174	(01)07290108240054(17)200501(10)9174(90)36098126(91)6070112	1-May-20
6070112	9271	(01)07290108240054(17)200501(10)9271(90)36098126(91)6070112	1-May-20
6070112	9272	(01)07290108240054(17)200501(10)9272(90)36098126(91)6070112	1-May-20
6070112	9273	(01)07290108240054(17)200501(10)9273(90)36098126(91)6070112	1-May-20
6070112	9274	(01)07290108240054(17)200501(10)9274(90)36098126(91)6070112	1-May-20
6070112	9275	(01)07290108240054(17)200601(10)9275(90)36098126(91)6070112	1-Jun-20
6070112	9276	(01)07290108240054(17)200601(10)9276(90)36098126(91)6070112	1-Jun-20
6070112	9277	(01)07290108240054(17)200601(10)9277(90)36098126(91)6070112	1-Jun-20
6070112	9278	(01)07290108240054(17)200601(10)9278(90)36098126(91)6070112	1-Jun-20



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	9323	(01)07290108240054(17)200601(10)9323(90)36098126(91)6070112	1-Jun-20
6070112	9324	(01)07290108240054(17)200601(10)9324(90)36098126(91)6070112	1-Jun-20
6070112	9325	(01)07290108240054(17)200601(10)9325(90)36098126(91)6070112	1-Jun-20
6070112	9326	(01)07290108240054(17)200601(10)9326(90)36098126(91)6070112	1-Jun-20
6070112	9327	(01)07290108240054(17)200601(10)9327(90)36098126(91)6070112	1-Jun-20
6070112	9328	(01)07290108240054(17)200601(10)9328(90)36098126(91)6070112	1-Jun-20
6070112	9329	(01)07290108240054(17)200601(10)9329(90)36098126(91)6070112	1-Jun-20
6070112	9330	(01)07290108240054(17)200601(10)9330(90)36098126(91)6070112	1-Jun-20
6070112	9394	(01)07290108240054(17)200501(10)9394	1-May-20
6070112	9500	(01)07290108240054(17)200701(10)9500(90)36098126(91)6070112	1-Jul-20
6070112	9501	(01)07290108240054(17)200701(10)9501(90)36098126(91)6070112	1-Jul-20
6070112	9502	(01)07290108240054(17)200701(10)9502(90)36098126(91)6070112	1-Jul-20
6070112	9503	(01)07290108240054(17)200701(10)9503(90)36098126(91)6070112	1-Jul-20
6070112	9504	(01)07290108240054(17)200701(10)9504(90)36098126(91)6070112	1-Jul-20
6070112	9505	(01)07290108240054(17)200701(10)9505(90)36098126(91)6070112	1-Jul-20
6070112	9506	(01)07290108240054(17)200701(10)9506(90)36098126(91)6070112	1-Jul-20
6070112	9507	(01)07290108240054(17)200701(10)9507(90)36098126(91)6070112	1-Jul-20
6070112	9508	(01)07290108240054(17)200801(10)9508(90)36098126(91)6070112	1-Aug-20
6070112	9509	(01)07290108240054(17)200801(10)9509(90)36098126(91)6070112	1-Aug-20
6070112	9510	(01)07290108240054(17)200801(10)9510(90)36098126(91)6070112	1-Aug-20
6070112	9511	(01)07290108240054(17)200801(10)9511(90)36098126(91)6070112	1-Aug-20
6070112	9512	(01)07290108240054(17)200801(10)9512(90)36098126(91)6070112	1-Aug-20
6070112	9513	(01)07290108240054(17)200801(10)9513(90)36098126(91)6070112	1-Aug-20
6070112	9514	(01)07290108240054(17)200801(10)9514(90)36098126(91)6070112	1-Aug-20
6070112	9515	(01)07290108240054(17)200801(10)9515(90)36098126(91)6070112	1-Aug-20
6070112	9529	(01)07290108240054(17)200801(10)9529(90)36098126(91)6070112	1-Aug-20
6070112	9530	(01)07290108240054(17)200801(10)9530(90)36098126(91)6070112	1-Aug-20
6070112	9531	(01)07290108240054(17)200801(10)9531(90)36098126(91)6070112	1-Aug-20
6070112	9532	(01)07290108240054(17)200801(10)9532(90)36098126(91)6070112	1-Aug-20
6070112	9533	(01)07290108240054(17)200801(10)9533(90)36098126(91)6070112	1-Aug-20
6070112	9534	(01)07290108240054(17)200801(10)9534(90)36098126(91)6070112	1-Aug-20
6070112	9535	(01)07290108240054(17)200801(10)9535(90)36098126(91)6070112	1-Aug-20
6070112	9537	(01)07290108240054(17)200801(10)9537(90)36098126(91)6070112	1-Aug-20
6070112	9560	(01)07290108240054(17)200901(10)9560(90)36098126(91)6070112	1-Sep-20
6070112	9561	(01)07290108240054(17)200901(10)9561(90)36098126(91)6070112	1-Sep-20
6070112	9562	(01)07290108240054(17)200901(10)9562(90)36098126(91)6070112	1-Sep-20
6070112	9563	(01)07290108240054(17)200901(10)9563(90)36098126(91)6070112	1-Sep-20



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	9564	(01)07290108240054(17)200901(10)9564(90)36098126(91)6070112	1-Sep-20
6070112	9565	(01)07290108240054(17)200901(10)9565(90)36098126(91)6070112	1-Sep-20
6070112	9566	(01)07290108240054(17)200901(10)9566(90)36098126(91)6070112	1-Sep-20
6070112	9567	(01)07290108240054(17)200901(10)9567(90)36098126(91)6070112	1-Sep-20
6070112	9568	(01)07290108240054(17)200901(10)9568(90)36098126(91)6070112	1-Sep-20
6070112	9569	(01)07290108240054(17)200901(10)9569(90)36098126(91)6070112	1-Sep-20
6070112	9570	(01)07290108240054(17)201001(10)9570(90)36098126(91)6070112	1-Oct-20
6070112	9571	(01)07290108240054(17)201001(10)9571(90)36098126(91)6070112	1-Oct-20
6070112	9572	(01)07290108240054(17)201001(10)9572(90)36098126(91)6070112	1-Oct-20
6070112	9657	(01)07290108240054(17)201001(10)9657(90)36098126(91)6070112	1-Oct-20
6070112	9658	(01)07290108240054(17)201001(10)9658(90)36098126(91)6070112	1-Oct-20
6070112	9659	(01)07290108240054(17)201001(10)9659(90)36098126(91)6070112	1-Oct-20
6070112	9660	(01)07290108240054(17)201001(10)9660(90)36098126(91)6070112	1-Oct-20
6070112	9661	(01)07290108240054(17)201101(10)9661(90)36098126(91)6070112	1-Nov-20
6070112	9662	(01)07290108240054(17)201101(10)9662(90)36098126(91)6070112	1-Nov-20
6070112	9663	(01)07290108240054(17)201101(10)9663(90)36098126(91)6070112	1-Nov-20
6070112	9664	(01)07290108240054(17)201101(10)9664(90)36098126(91)6070112	1-Nov-20
6070112	9665	(01)07290108240054(17)201101(10)9665(90)36098126(91)6070112	1-Nov-20
6070112	9666	(01)07290108240054(17)201101(10)9666(90)36098126(91)6070112	1-Nov-20
6070112	9667	(01)07290108240054(17)201101(10)9667(90)36098126(91)6070112	1-Nov-20
6070112	9668	(01)07290108240054(17)201101(10)9668(90)36098126(91)6070112	1-Nov-20
6070112	9717	(01)07290108240054(17)201101(10)9717(90)36098126(91)6070112	1-Nov-20
6070112	9718	(01)07290108240054(17)201101(10)9718(90)36098126(91)6070112	1-Nov-20
6070112	9719	(01)07290108240054(17)201101(10)9719(90)36098126(91)6070112	1-Nov-20
6070112	9720	(01)07290108240054(17)201101(10)9720(90)36098126(91)6070112	1-Nov-20
6070112	9721	(01)07290108240054(17)201101(10)9721(90)36098126(91)6070112	1-Nov-20
6070112	9722	(01)07290108240054(17)201201(10)9722(90)36098126(91)6070112	1-Dec-20
6070112	9723	(01)07290108240054(17)201201(10)9723(90)36098126(91)6070112	1-Dec-20
6070112	9724	(01)07290108240054(17)201201(10)9724(90)36098126(91)6070112	1-Dec-20
6070112	9725	(01)07290108240054(17)201201(10)9725(90)36098126(91)6070112	1-Dec-20
6070112	9726	(01)07290108240054(17)201201(10)9726(90)36098126(91)6070112	1-Dec-20
6070112	9727	(01)07290108240054(17)201201(10)9727(90)36098126(91)6070112	1-Dec-20
6070112	9728	(01)07290108240054(17)201201(10)9728(90)36098126(91)6070112	1-Dec-20
6070112	9729	(01)07290108240054(17)201201(10)9729(90)36098126(91)6070112	1-Dec-20
6070112	9730	(01)07290108240054(17)201201(10)9730(90)36098126(91)6070112	1-Dec-20
6070112	9810	(01)07290108240054(17)200901(10)9810(90)36098126(91)6070112	1-Sep-20
6070112	9841	(01)07290108240054(17)201201(10)9841(90)36098126(91)6070112	1-Dec-20



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	9842	(01)07290108240054(17)201201(10)9842(90)36098126(91)6070112	1-Dec-20
6070112	9843	(01)07290108240054(17)201201(10)9843(90)36098126(91)6070112	1-Dec-20
6070112	9844	(01)07290108240054(17)201201(10)9844(90)36098126(91)6070112	1-Dec-20
6070112	9845	(01)07290108240054(17)201201(10)9845(90)36098126(91)6070112	1-Dec-20
6070112	9846	(01)07290108240054(17)201201(10)9846(90)36098126(91)6070112	1-Dec-20
6070112	9847	(01)07290108240054(17)201201(10)9847(90)36098126(91)6070112	1-Dec-20
6070112	9848	(01)07290108240054(17)201201(10)9848(90)36098126(91)6070112	1-Dec-20
6070112	9849	(01)07290108240054(17)201201(10)9849(90)36098126(91)6070112	1-Dec-20
6070112	9964	(01)07290108240054(17)210101(10)9964(90)36098126(91)6070112	1-Jan-21
6070112	9965	(01)07290108240054(17)210101(10)9965(90)36098126(91)6070112	1-Jan-21
6070112	9966	(01)07290108240054(17)210101(10)9966(90)36098126(91)6070112	1-Jan-21
6070112	9967	(01)07290108240054(17)210101(10)9967(90)36098126(91)6070112	1-Jan-21
6070112	9968	(01)07290108240054(17)210101(10)9968(90)36098126(91)6070112	1-Jan-21
6070112	9969	(01)07290108240054(17)210101(10)9969(90)36098126(91)6070112	1-Jan-21
6070112	9970	(01)07290108240054(17)210101(10)9970(90)36098126(91)6070112	1-Jan-21
6070112	9971	(01)07290108240054(17)210101(10)9971(90)36098126(91)6070112	1-Jan-21
6070112	9972	(01)07290108240054(17)210101(10)9972(90)36098126(91)6070112	1-Jan-21
6070112	9973	(01)07290108240054(17)210101(10)9973(90)36098126(91)6070112	1-Jan-21
6070112	9974	(01)07290108240054(17)210101(10)9974(90)36098126(91)6070112	1-Jan-21
6070112	9975	(01)07290108240054(17)210101(10)9975(90)36098126(91)6070112	1-Jan-21
6070112	9976	(01)07290108240054(17)210101(10)9976(90)36098126(91)6070112	1-Jan-21
6070112	9978	(01)07290108240054(17)210101(10)9978(90)36098126(91)6070112	1-Jan-21
6070112	9979	(01)07290108240054(17)210101(10)9979(90)36098126(91)6070112	1-Jan-21
6070112	9980	(01)07290108240054(17)210101(10)9980(90)36098126(91)6070112	1-Jan-21
6070112	A001	(01)07290108240054(17)210201(10)A001(90)36098126(91)6070112	1-Feb-21
6070112	A002	(01)07290108240054(17)210201(10)A002(90)36098126(91)6070112	1-Feb-21
6070112	A003	(01)07290108240054(17)210201(10)A003(90)36098126(91)6070112	1-Feb-21
6070112	A004	(01)07290108240054(17)210201(10)A004(90)36098126(91)6070112	1-Feb-21
6070112	A005	(01)07290108240054(17)210201(10)A005(90)36098126(91)6070112	1-Feb-21
6070112	A006	(01)07290108240054(17)210301(10)A006(90)36098126(91)6070112	1-Mar-21
6070112	A007	(01)07290108240054(17)210301(10)A007(90)36098126(91)6070112	1-Mar-21
6070112	A008	(01)07290108240054(17)210301(10)A008(90)36098126(91)6070112	1-Mar-21
6070112	A009	(01)07290108240054(17)210301(10)A009(90)36098126(91)6070112	1-Mar-21
6070112	A010	(01)07290108240054(17)210301(10)A010(90)36098126(91)6070112	1-Mar-21
6070112	A011	(01)07290108240054(17)210301(10)A011(90)36098126(91)6070112	1-Mar-21
6070112	A012	(01)07290108240054(17)210301(10)A012(90)36098126(91)6070112	1-Mar-21
6070112	A013	(01)07290108240054(17)210301(10)A013(90)36098126(91)6070112	1-Mar-21



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	A014	(01)07290108240054(17)210301(10)A014(90)36098126(91)6070112	1-Mar-21
6070112	A015	(01)07290108240054(17)210301(10)A015(90)36098126(91)6070112	1-Mar-21
6070112	A016	(01)07290108240054(17)210301(10)A016(90)36098126(91)6070112	1-Mar-21
6070112	A017	(01)07290108240054(17)210301(10)A017(90)36098126(91)6070112	1-Mar-21
6070112	A018	(01)07290108240054(17)210101(10)A018(90)36098126(91)6070112	1-Jan-21
6070112	A019	(01)07290108240054(17)210101(10)A019(90)36098126(91)6070112	1-Jan-21
6070112	A020	(01)07290108240054(17)210201(10)A020(90)36098126(91)6070112	1-Feb-21
6070112	A021	(01)07290108240054(17)210201(10)A021(90)36098126(91)6070112	1-Feb-21
6070112	A022	(01)07290108240054(17)210201(10)A022(90)36098126(91)6070112	1-Feb-21
6070112	A023	(01)07290108240054(17)210201(10)A023(90)36098126(91)6070112	1-Feb-21
6070112	A024	(01)07290108240054(17)210201(10)A024(90)36098126(91)6070112	1-Feb-21
6070112	A025	(01)07290108240054(17)210201(10)A025(90)36098126(91)6070112	1-Feb-21
6070112	A026	(01)07290108240054(17)210201(10)A026(90)36098126(91)6070112	1-Feb-21
6070112	A027	(01)07290108240054(17)210201(10)A027(90)36098126(91)6070112	1-Feb-21
6070112	A028	(01)07290108240054(17)210201(10)A028(90)36098126(91)6070112	1-Feb-21
6070112	A029	(01)07290108240054(17)210201(10)A029(90)36098126(91)6070112	1-Feb-21
6070112	A030	(01)07290108240054(17)210201(10)A030(90)36098126(91)6070112	1-Feb-21
6070112	A031	(01)07290108240054(17)210201(10)A031(90)36098126(91)6070112	1-Feb-21
6070112	A032	(01)07290108240054(17)210201(10)A032(90)36098126(91)6070112	1-Feb-21
6070112	A033	(01)07290108240054(17)210201(10)A033(90)36098126(91)6070112	1-Feb-21
6070112	A034	(01)07290108240054(17)210201(10)A034(90)36098126(91)6070112	1-Feb-21
6070112	A106	(01)07290108240054(17)210301(10)A106(90)36098126(91)6070112	1-Mar-21
6070112	A107	(01)07290108240054(17)210301(10)A107(90)36098126(91)6070112	1-Mar-21
6070112	A108	(01)07290108240054(17)210301(10)A108(90)36098126(91)6070112	1-Mar-21
6070112	A109	(01)07290108240054(17)210301(10)A109(90)36098126(91)6070112	1-Mar-21
6070112	A110	(01)07290108240054(17)210301(10)A110(90)36098126(91)6070112	1-Mar-21
6070112	A111	(01)07290108240054(17)210301(10)A111(90)36098126(91)6070112	1-Mar-21
6070112	A112	(01)07290108240054(17)210401(10)A112(90)36098126(91)6070112	1-Apr-21
6070112	A113	(01)07290108240054(17)210401(10)A113(90)36098126(91)6070112	1-Apr-21
6070112	A114	(01)07290108240054(17)210401(10)A114(90)36098126(91)6070112	1-Apr-21
6070112	A115	(01)07290108240054(17)210401(10)A115(90)36098126(91)6070112	1-Apr-21
6070112	A116	(01)07290108240054(17)210401(10)A116(90)36098126(91)6070112	1-Apr-21
6070112	A117	(01)07290108240054(17)210401(10)A117(90)36098126(91)6070112	1-Apr-21
6070112	A118	(01)07290108240054(17)210401(10)A118(90)36098126(91)6070112	1-Apr-21
6070112	A119	(01)07290108240054(17)210401(10)A119(90)36098126(91)6070112	1-Apr-21
6070112	A120	(01)07290108240054(17)210401(10)A120(90)36098126(91)6070112	1-Apr-21
6070112	A121	(01)07290108240054(17)210401(10)A121(90)36098126(91)6070112	1-Apr-21



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	A122	(01)07290108240054(17)210401(10)A122(90)36098126(91)6070112	1-Apr-21
6070112	A254	(01)07290108240054(17)210201(10)A254	1-Feb-21
6070112	A258	(01)07290108240054(17)210301(10)A258(90)36098126(91)6070112	1-Mar-21
6070112	A321	(01)07290108240054(17)210401(10)A321(90)36098126(91)6070112	1-Apr-21
6070112	A322	(01)07290108240054(17)210401(10)A322(90)36098126(91)6070112	1-Apr-21
6070112	A323	(01)07290108240054(17)210401(10)A323(90)36098126(91)6070112	1-Apr-21
6070112	A324	(01)07290108240054(17)210501(10)A324(90)36098126(91)6070112	1-May-21
6070112	A325	(01)07290108240054(17)210501(10)A325(90)36098126(91)6070112	1-May-21
6070112	A326	(01)07290108240054(17)210501(10)A326(90)36098126(91)6070112	1-May-21
6070112	A327	(01)07290108240054(17)210501(10)A327(90)36098126(91)6070112	1-May-21
6070112	A328	(01)07290108240054(17)210501(10)A328(90)36098126(91)6070112	1-May-21
6070112	A329	(01)07290108240054(17)210501(10)A329(90)36098126(91)6070112	1-May-21
6070112	A330	(01)07290108240054(17)210501(10)A330(90)36098126(91)6070112	1-May-21
6070112	A331	(01)07290108240054(17)210501(10)A331(90)36098126(91)6070112	1-May-21
6070112	A332	(01)07290108240054(17)210501(10)A332(90)36098126(91)6070112	1-May-21
6070112	A333	(01)07290108240054(17)210501(10)A333(90)36098126(91)6070112	1-May-21
6070112	A334	(01)07290108240054(17)210501(10)A334(90)36098126(91)6070112	1-May-21
6070112	A335	(01)07290108240054(17)210501(10)A335(90)36098126(91)6070112	1-May-21
6070112	A336	(01)07290108240054(17)210501(10)A336(90)36098126(91)6070112	1-May-21
6070112	A337	(01)07290108240054(17)210501(10)A337(90)36098126(91)6070112	1-May-21
6070112	A338	(01)07290108240054(17)210601(10)A338(90)36098126(91)6070112	1-Jun-21
6070112	A339	(01)07290108240054(17)210601(10)A339(90)36098126(91)6070112	1-Jun-21
6070112	A340	(01)07290108240054(17)210601(10)A340(90)36098126(91)6070112	1-Jun-21
6070112	A343	(01)07290108240054(17)210401(10)A343(90)36098126(91)6070112	1-Apr-21
6070112	A344	(01)07290108240054(17)210401(10)A344(90)36098126(91)6070112	1-Apr-21
6070112	A345	(01)07290108240054(17)210401(10)A345(90)36098126(91)6070112	1-Apr-21
6070112	A346	(01)07290108240054(17)210401(10)A346(90)36098126(91)6070112	1-Apr-21
6070112	A347	(01)07290108240054(17)210401(10)A347(90)36098126(91)6070112	1-Apr-21
6070112	A348	(01)07290108240054(17)210401(10)A348(90)36098126(91)6070112	1-Apr-21
6070112	A349	(01)07290108240054(17)210401(10)A349(90)36098126(91)6070112	1-Apr-21
6070112	A350	(01)07290108240054(17)210401(10)A350(90)36098126(91)6070112	1-Apr-21
6070112	A351	(01)07290108240054(17)210401(10)A351(90)36098126(91)6070112	1-Apr-21
6070112	A352	(01)07290108240054(17)210401(10)A352(90)36098126(91)6070112	1-Apr-21
6070112	A353	(01)07290108240054(17)210401(10)A353(90)36098126(91)6070112	1-Apr-21
6070112	A354	(01)07290108240054(17)210401(10)A354(90)36098126(91)6070112	1-Apr-21
6070112	A355	(01)07290108240054(17)210401(10)A355(90)36098126(91)6070112	1-Apr-21
6070112	A356	(01)07290108240054(17)210501(10)A356(90)36098126(91)6070112	1-May-21



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	A357	(01)07290108240054(17)210501(10)A357(90)36098126(91)6070112	1-May-21
6070112	A358	(01)07290108240054(17)210501(10)A358(90)36098126(91)6070112	1-May-21
6070112	A359	(01)07290108240054(17)210501(10)A359(90)36098126(91)6070112	1-May-21
6070112	A360	(01)07290108240054(17)210501(10)A360(90)36098126(91)6070112	1-May-21
6070112	A361	(01)07290108240054(17)210501(10)A361(90)36098126(91)6070112	1-May-21
6070112	A362	(01)07290108240054(17)210501(10)A362(90)36098126(91)6070112	1-May-21
6070112	A519	(01)07290108240054(17)210601(10)A519(90)36098126(91)6070112	1-Jun-21
6070112	A520	(01)07290108240054(17)210601(10)A520(90)36098126(91)6070112	1-Jun-21
6070112	A521	(01)07290108240054(17)210601(10)A521(90)36098126(91)6070112	1-Jun-21
6070112	A522	(01)07290108240054(17)210601(10)A522(90)36098126(91)6070112	1-Jun-21
6070112	A523	(01)07290108240054(17)210601(10)A523(90)36098126(91)6070112	1-Jun-21
6070112	A524	(01)07290108240054(17)210601(10)A524(90)36098126(91)6070112	1-Jun-21
6070112	A525	(01)07290108240054(17)210601(10)A525(90)36098126(91)6070112	1-Jun-21
6070112	A526	(01)07290108240054(17)210601(10)A526(90)36098126(91)6070112	1-Jun-21
6070112	A527	(01)07290108240054(17)210601(10)A527(90)36098126(91)6070112	1-Jun-21
6070112	A528	(01)07290108240054(17)210601(10)A528(90)36098126(91)6070112	1-Jun-21
6070112	A529	(01)07290108240054(17)210601(10)A529(90)36098126(91)6070112	1-Jun-21
6070112	A530	(01)07290108240054(17)210601(10)A530(90)36098126(91)6070112	1-Jun-21
6070112	A531	(01)07290108240054(17)210601(10)A531(90)36098126(91)6070112	1-Jun-21
6070112	A532	(01)07290108240054(17)210601(10)A532(90)36098126(91)6070112	1-Jun-21
6070112	A533	(01)07290108240054(17)210601(10)A533(90)36098126(91)6070112	1-Jun-21
6070112	A534	(01)07290108240054(17)210601(10)A534(90)36098126(91)6070112	1-Jun-21
6070112	A535	(01)07290108240054(17)210601(10)A535(90)36098126(91)6070112	1-Jun-21
6070112	A536	(01)07290108240054(17)210701(10)A536(90)36098126(91)6070112	1-Jul-21
6070112	A537	(01)07290108240054(17)210701(10)A537(90)36098126(91)6070112	1-Jul-21
6070112	A538	(01)07290108240054(17)210701(10)A538(90)36098126(91)6070112	1-Jul-21
6070112	A539	(01)07290108240054(17)210701(10)A539(90)36098126(91)6070112	1-Jul-21
6070112	A540	(01)07290108240054(17)210701(10)A540(90)36098126(91)6070112	1-Jul-21
6070112	A541	(01)07290108240054(17)210701(10)A541(90)36098126(91)6070112	1-Jul-21
6070112	A542	(01)07290108240054(17)210701(10)A542(90)36098126(91)6070112	1-Jul-21
6070112	A543	(01)07290108240054(17)210701(10)A543(90)36098126(91)6070112	1-Jul-21
6070112	A544	(01)07290108240054(17)210701(10)A544(90)36098126(91)6070112	1-Jul-21
6070112	A545	(01)07290108240054(17)210701(10)A545(90)36098126(91)6070112	1-Jul-21
6070112	A546	(01)07290108240054(17)210701(10)A546(90)36098126(91)6070112	1-Jul-21
6070112	A547	(01)07290108240054(17)210701(10)A547(90)36098126(91)6070112	1-Jul-21
6070112	A548	(01)07290108240054(17)210701(10)A548(90)36098126(91)6070112	1-Jul-21
6070112	A549	(01)07290108240054(17)210701(10)A549(90)36098126(91)6070112	1-Jul-21



West Pharmaceutical Services, Inc.
 530 Herman O. West Drive
 Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	A550	(01)07290108240054(17)210701(10)A550(90)36098126(91)6070112	1-Jul-21
6070112	A551	(01)07290108240054(17)210701(10)A551(90)36098126(91)6070112	1-Jul-21
6070112	A552	(01)07290108240054(17)210701(10)A552(90)36098126(91)6070112	1-Jul-21
6070112	A553	(01)07290108240054(17)210701(10)A553(90)36098126(91)6070112	1-Jul-21
6070112	A554	(01)07290108240054(17)210701(10)A554(90)36098126(91)6070112	1-Jul-21
6070112	A555	(01)07290108240054(17)210701(10)A555(90)36098126(91)6070112	1-Jul-21
6070112	A556	(01)07290108240054(17)210701(10)A556(90)36098126(91)6070112	1-Jul-21
6070112	A557	(01)07290108240054(17)210701(10)A557(90)36098126(91)6070112	1-Jul-21
6070112	A558	(01)07290108240054(17)210701(10)A558(90)36098126(91)6070112	1-Jul-21
6070112	A559	(01)07290108240054(17)210701(10)A559(90)36098126(91)6070112	1-Jul-21
6070112	A560	(01)07290108240054(17)210701(10)A560(90)36098126(91)6070112	1-Jul-21
6070112	A561	(01)07290108240054(17)210701(10)A561(90)36098126(91)6070112	1-Jul-21
6070112	A562	(01)07290108240054(17)210701(10)A562(90)36098126(91)6070112	1-Jul-21
6070112	A563	(01)07290108240054(17)210701(10)A563(90)36098126(91)6070112	1-Jul-21
6070112	A564	(01)07290108240054(17)210701(10)A564(90)36098126(91)6070112	1-Jul-21
6070112	A565	(01)07290108240054(17)210701(10)A565(90)36098126(91)6070112	1-Jul-21
6070112	A566	(01)07290108240054(17)210701(10)A566(90)36098126(91)6070112	1-Jul-21
6070112	A567	(01)07290108240054(17)210801(10)A567(90)36098126(91)6070112	1-Aug-21
6070112	A568	(01)07290108240054(17)210701(10)A568(90)36098126(91)6070112	1-Jul-21
6070112	A569	(01)07290108240054(17)210701(10)A569(90)36098126(91)6070112	1-Jul-21
6070112	A570	(01)07290108240054(17)210701(10)A570(90)36098126(91)6070112	1-Jul-21
6070112	A571	(01)07290108240054(17)210701(10)A571(90)36098126(91)6070112	1-Jul-21
6070112	A572	(01)07290108240054(17)210701(10)A572(90)36098126(91)6070112	1-Jul-21
6070112	A573	(01)07290108240054(17)210801(10)A573(90)36098126(91)6070112	1-Aug-21
6070112	A574	(01)07290108240054(17)210801(10)A574(90)36098126(91)6070112	1-Aug-21
6070112	A575	(01)07290108240054(17)210801(10)A575(90)36098126(91)6070112	1-Aug-21
6070112	A576	(01)07290108240054(17)210801(10)A576(90)36098126(91)6070112	1-Aug-21
6070112	A577	(01)07290108240054(17)210801(10)A577(90)36098126(91)6070112	1-Aug-21
6070112	A578	(01)07290108240054(17)210801(10)A578(90)36098126(91)6070112	1-Aug-21
6070112	A579	(01)07290108240054(17)210801(10)A579(90)36098126(91)6070112	1-Aug-21
6070112	A580	(01)07290108240054(17)210801(10)A580(90)36098126(91)6070112	1-Aug-21
6070112	A581	(01)07290108240054(17)210801(10)A581(90)36098126(91)6070112	1-Aug-21
6070112	A582	(01)07290108240054(17)210801(10)A582(90)36098126(91)6070112	1-Aug-21
6070112	A630	(01)07290108240054(17)210701(10)A630	1-Jul-21
6070112	A813	(01)07290108240054(17)210901(10)A813(90)36098126(91)6070112	1-Sep-21
6070112	A814	(01)07290108240054(17)210901(10)A814(90)36098126(91)6070112	1-Sep-21
6070112	A815	(01)07290108240054(17)210901(10)A815(90)36098126(91)6070112	1-Sep-21



West Pharmaceutical Services, Inc.
 530 Herman O. West Drive
 Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	A816	(01)07290108240054(17)210901(10)A816(90)36098126(91)6070112	1-Sep-21
6070112	A817	(01)07290108240054(17)210801(10)A817(90)36098126(91)6070112	1-Aug-21
6070112	A818	(01)07290108240054(17)210901(10)A818(90)36098126(91)6070112	1-Sep-21
6070112	A819	(01)07290108240054(17)210801(10)A819(90)36098126(91)6070112	1-Aug-21
6070112	A820	(01)07290108240054(17)210901(10)A820(90)36098126(91)6070112	1-Sep-21
6070112	A821	(01)07290108240054(17)210901(10)A821(90)36098126(91)6070112	1-Sep-21
6070112	A822	(01)07290108240054(17)210901(10)A822(90)36098126(91)6070112	1-Sep-21
6070112	A823	(01)07290108240054(17)210901(10)A823(90)36098126(91)6070112	1-Sep-21
6070112	A824	(01)07290108240054(17)210901(10)A824(90)36098126(91)6070112	1-Sep-21
6070112	A825	(01)07290108240054(17)210901(10)A825(90)36098126(91)6070112	1-Sep-21
6070112	A826	(01)07290108240054(17)210901(10)A826(90)36098126(91)6070112	1-Sep-21
6070112	A840	(01)07290108240054(17)210901(10)A840(90)36098126(91)6070112	1-Sep-21
6070112	A841	(01)07290108240054(17)210901(10)A841(90)36098126(91)6070112	1-Sep-21
6070112	A842	(01)07290108240054(17)210901(10)A842(90)36098126(91)6070112	1-Sep-21
6070112	A843	(01)07290108240054(17)210901(10)A843(90)36098126(91)6070112	1-Sep-21
6070112	A844	(01)07290108240054(17)210901(10)A844(90)36098126(91)6070112	1-Sep-21
6070112	A845	(01)07290108240054(17)210901(10)A845(90)36098126(91)6070112	1-Sep-21
6070112	A898	(01)07290108240054(17)211001(10)A898(90)36098126(91)6070112	1-Oct-21
6070112	A899	(01)07290108240054(17)211001(10)A899(90)36098126(91)6070112	1-Oct-21
6070112	A900	(01)07290108240054(17)211001(10)A900(90)36098126(91)6070112	1-Oct-21
6070112	A901	(01)07290108240054(17)211001(10)A901(90)36098126(91)6070112	1-Oct-21
6070112	A902	(01)07290108240054(17)211001(10)A902(90)36098126(91)6070112	1-Oct-21
6070112	A903	(01)07290108240054(17)211001(10)A903(90)36098126(91)6070112	1-Oct-21
6070112	A904	(01)07290108240054(17)211001(10)A904(90)36098126(91)6070112	1-Oct-21
6070112	A905	(01)07290108240054(17)211001(10)A905(90)36098126(91)6070112	1-Oct-21
6070112	A906	(01)07290108240054(17)211001(10)A906(90)36098126(91)6070112	1-Oct-21
6070112	A907	(01)07290108240054(17)211001(10)A907(90)36098126(91)6070112	1-Oct-21
6070112	A908	(01)07290108240054(17)211001(10)A908(90)36098126(91)6070112	1-Oct-21
6070112	A909	(01)07290108240054(17)211001(10)A909(90)36098126(91)6070112	1-Oct-21
6070112	A910	(01)07290108240054(17)211001(10)A910(90)36098126(91)6070112	1-Oct-21
6070112	A911	(01)07290108240054(17)211001(10)A911(90)36098126(91)6070112	1-Oct-21
6070112	A915	(01)07290108240054(17)211001(10)A915(90)36098126(91)6070112	1-Oct-21
6070112	A916	(01)07290108240054(17)211001(10)A916(90)36098126(91)6070112	1-Oct-21
6070112	A917	(01)07290108240054(17)211001(10)A917(90)36098126(91)6070112	1-Oct-21
6070112	A918	(01)07290108240054(17)211001(10)A918(90)36098126(91)6070112	1-Oct-21
6070112	A919	(01)07290108240054(17)211001(10)A919(90)36098126(91)6070112	1-Oct-21
6070112	A920	(01)07290108240054(17)211001(10)A920(90)36098126(91)6070112	1-Oct-21



West Pharmaceutical Services, Inc.
530 Herman O. West Drive
Exton, PA 19341
www.westpharma.com

Vial2Bag DC 20mm			
Part Number	Lot	UDI	Expiry Date
6070112	A921	(01)07290108240054(17)211001(10)A921(90)36098126(91)6070112	1-Oct-21
6070112	A982	(01)07290108240054(17)211001(10)A982(90)36098126(91)6070112	1-Oct-21
6070112	A983	(01)07290108240054(17)211001(10)A983(90)36098126(91)6070112	1-Oct-21
6070112	A984	(01)07290108240054(17)211001(10)A984(90)36098126(91)6070112	1-Oct-21
6070112	A985	(01)07290108240054(17)211101(10)A985(90)36098126(91)6070112	1-Nov-21
6070112	A986	(01)07290108240054(17)211101(10)A986(90)36098126(91)6070112	1-Nov-21
6070112	A987	(01)07290108240054(17)211101(10)A987(90)36098126(91)6070112	1-Nov-21
6070112	A988	(01)07290108240054(17)211101(10)A988(90)36098126(91)6070112	1-Nov-21
6070112	C015	(01)07290108240054(17)211101(10)C015(90)36098126(91)6070112	1-Nov-21
6070112	C016	(01)07290108240054(17)211101(10)C016(90)36098126(91)6070112	1-Nov-21
6070112	C017	(01)07290108240054(17)211101(10)C017(90)36098126(91)6070112	1-Nov-21
6070112	C033	(01)07290108240054(17)211101(10)C033(90)36098126(91)6070112	1-Nov-21
6070112	C034	(01)07290108240054(17)211101(10)C034(90)36098126(91)6070112	1-Nov-21
6070112	C035	(01)07290108240054(17)211101(10)C035(90)36098126(91)6070112	1-Nov-21
6070112	C036	(01)07290108240054(17)211101(10)C036(90)36098126(91)6070112	1-Nov-21
6070112	C037	(01)07290108240054(17)211101(10)C037(90)36098126(91)6070112	1-Nov-21
6070112	C038	(01)07290108240054(17)211101(10)C038(90)36098126(91)6070112	1-Nov-21
6070112	C039	(01)07290108240054(17)211101(10)C039(90)36098126(91)6070112	1-Nov-21
6070112	C040	(01)07290108240054(17)211101(10)C040(90)36098126(91)6070112	1-Nov-21
6070112	C041	(01)07290108240054(17)211101(10)C041(90)36098126(91)6070112	1-Nov-21
6070112	C042	(01)07290108240054(17)211101(10)C042(90)36098126(91)6070112	1-Nov-21
6070112	C043	(01)07290108240054(17)211101(10)C043(90)36098126(91)6070112	1-Nov-21
6070112	C044	(01)07290108240054(17)211101(10)C044(90)36098126(91)6070112	1-Nov-21
6070112	C045	(01)07290108240054(17)211101(10)C045(90)36098126(91)6070112	1-Nov-21
6070112	C046	(01)07290108240054(17)211101(10)C046(90)36098126(91)6070112	1-Nov-21
6070112	C047	(01)07290108240054(17)211101(10)C047(90)36098126(91)6070112	1-Nov-21
6070112	C048	(01)07290108240054(17)211101(10)C048(90)36098126(91)6070112	1-Nov-21
6070112	C049	(01)07290108240054(17)211101(10)C049(90)36098126(91)6070112	1-Nov-21
6070112	C050	(01)07290108240054(17)211101(10)C050(90)36098126(91)6070112	1-Nov-21
6070112	C051	(01)07290108240054(17)211101(10)C051(90)36098126(91)6070112	1-Nov-21
6070112	C052	(01)07290108240054(17)211101(10)C052(90)36098126(91)6070112	1-Nov-21
6070112	C053	(01)07290108240054(17)211101(10)C053(90)36098126(91)6070112	1-Nov-21
6070112	C054	(01)07290108240054(17)211101(10)C054(90)36098126(91)6070112	1-Nov-21
6070112	C055	(01)07290108240054(17)211101(10)C055(90)36098126(91)6070112	1-Nov-21
6070112	C056	(01)07290108240054(17)211101(10)C056(90)36098126(91)6070112	1-Nov-21
6070112C	8344	(01)07290108241006(17)190501(10)8344(90)36098138(91)6070112C	1-May-19

NOTE: The highlighted items indicate newly added affected lots