

4

3

2

1

PRIMARY BAG –  
STERILIZABLEBAG – Heat Sealed

SECONDARY BAG –  
STERILIZABLEBAG – Heat Sealed

TERTIARY CLEAR  
CR100 PE BAG – Heat Sealed

TERTIARY FOIL MOISTURE  
BARRIER BAG – Heat Sealed

SERUM

LYO

PLASTIC CARTON  
15 1/4" X 12 1/4" X 2 1/4" O.D.  
(1) DOUBLE STERILIZABLEBAG/CARTON

LABEL

LABEL


LABEL /  
USE BY DATE

STERILE  
CLAIM

LABEL /  
USE BY DATE

STERILE  
CLAIM

CARTON LABEL



REV.	ZONE	DATE	DESCRIPTION	DR.	APPRD.
0	GEN	25 JAN 2018	PR 53055	HG	MM
REVISIONS					

TOLERANCES UNLESS NOTED DECIMAL METRIC .XX ± .01    .X ± .3 .XXX ± .003    .XX ± .08 .XXXX ± .0015    .XXX ± .038		APPROVALS DWN: H. Giang CHKD: M. Habash DWG. APPVD: M. Mariani	DATE 17 JAN 2018 24 JAN 2018 25 JAN 2018		TITLE PACKAGE ASSEMBLY ILLUSTRATION ELASTOMER RU SMALL CARTON
FRACTIONAL INCH: ± 1/64 ANGULAR: ± 1° SURFACE FINISH: ✓ THE INFORMATION ON THIS DRAWING IS PROPRIETARY. THIS INFORMATION SHALL NOT BE DISCLOSED, REPRODUCED IN ANY FORM OR USED FOR ANY PURPOSE WITHOUT PRIOR WRITTEN CONSENT OF WEST PHARMACEUTICAL SERVICES. DO NOT SCALE DRAWING	TOOL ENGR: N/A PROD. ENGR: N/A APPRD: N/A SAP Number: N/A	DATE: N/A DATE: N/A DATE: N/A	MATL: N/A DRAWING AND REV. NO: 1841-022 Rev 0 SCALE: N/A    WEIGHT: N/A    SHEET 1 OF 1		